

Singfiller®

sterile injection
safety and natural
strong resistance


CE
1023




Cross-linked sodium Hyaluronate for Plastic Surgery

APPLICATION

- 1 facial beauty injection
- 2 improve skin appearance
- 3 structure, reduce wrinkle, prevent skin aging
- 4 improve labial shape
- 5 facial tissue augmentation



- ① Nasolabial folds
- ② Frown lines
- ③ Wrinkles across the forehead
- ④ Neck lines
- ⑤ Nose
- ⑥ Chin
- ⑦ Lip
- ⑧ Earlobe
- ⑨ Temple
- ⑩ Cheek




Deep	Deeper	Sub-Q
20mg/ml(HA) 0.28~0.50mm	20mg/ml(HA) 0.50~1.00mm	20mg/ml(HA) 1.25~2.00 mm
Deep facial wrinkles and folds ①⑥⑦⑧⑨⑩ deep layer of dermis	Deep facial wrinkle and folds ①⑨⑩ deep dermis or subcutaneous tissue shallow	the local depression, soft tissue filler. And breast enhance , Breast butt deep dermis or subcutaneous tissue shallow
1mL (CE) / 2mL (CE)	1mL (CE) / 2mL (CE) / 10mL(No CE)	2mL(No CE) / 20mL(No CE)
0~30°C (32-86°F)	0~30°C (32-86°F)	0~30°C (32-86°F)
Single	Single	Single
6 - 9 months	6-12months	6-12months
26G*1 27G*2	26G*1 26G*2	/ /


NASOLABIAL INJECTION



MECHANISM OF ACTION



Singfiller® shows equivalent results on Suitability, Safety & Effectivity.



Data from: "Summary of Clinical investigation for Cross-linked Sodium Hyaluronate Gel for Plastic Surgery"



Singfiller®

Cross-linked sodium Hyaluronate for Plastic Surgery



www.singclean.net
www.singfiller.org

Hangzhou Singclean Medical Products Co., Ltd
Add: No. 125 (E), 10th Street, Hangzhou Economic and
Technological Development Zone, ZheJiang, China
P.C.: 310018
Tel: +86-571-63431868 Fax: +86-571-63431886
E-mail: sales@hzxhe.com
Web: www.singclean.net www.singfiller.org